

SM1007

Thin Cylinder

Bench-mounted machine to allow students to do stress and strain tests on a thin-walled cylinder

Works with
VDAS®



- Ideal for student use and classroom demonstrations
- For comprehensive analysis of the stresses and strains in a thin-walled cylinder, under internal pressure
- Introduces stress and strain, and how to measure and analyse them
- Includes experiments to find Young's modulus and Poisson's ratio
- Closed-end and open-end conditions to allow circumferential or biaxial stress tests
- High-quality electrical-resistance precision strain gauges measure cylinder strains
- Includes built-in microprocessor-controlled display of strain measurements
- Mounted on a sturdy base, to form a compact product ideal for use on a workbench
- Self-contained, hand-operated hydraulic pressurising system for accurate pressure control
- Completely self-contained – needs no other parts

- TecEquipment Ltd, Bonsall Street, Long Eaton, Nottingham NG10 2AN, UK
- T +44 115 972 2611 • F +44 115 973 1520 • E info@tecequipment.com • W www.tecequipment.com
- An ISO 9001 certified company
- VDAS is a registered trademark of TecEquipment Ltd

SM1007

Thin Cylinder

Description

TecEquipment's Thin Cylinder apparatus allows students to perform experiments that examine stress and strain in a thin-walled cylinder. It clearly shows the principles, theories and analytical techniques, and provides effective, practical support to studies.

A sturdy base contains all parts of the Thin Cylinder apparatus. This forms a compact product, ideal for use on a workbench.

The apparatus consists of a thin-walled aluminium cylinder, held in a robust frame. The frame holds the cylinder so that it is free to move along its axis. The cylinder contains oil. To stress the cylinder, students use the hydraulic hand-pump to pressurise the oil. Strain gauges on the cylinder surface measure strain, while a gauge and electronic sensor measure hydraulic pressure.

Students can measure strains with the cylinder in two 'end conditions':

- Open end: the cylinder has no axial load, so there is no direct axial stress.

or

- Closed end: the cylinder has axial loads, so there is direct axial stress.

Students use a hand-wheel on the frame to set these end conditions.

To perform experiments, students choose either closed or open-end conditions. They set the gauges to zero and use the pump to pressurise the cylinder. They take readings at several stages while they increase the pressure. The results can be taken by hand using the in-built display and pressure gauge, and results plotted by hand. Alternatively, students can use TecEquipment's optional Versatile Data Acquisition System (VDAS®) to capture the data and plot the relevant graphs and export data. They then compare their results with calculations made using stress and strain theory.

A user guide is supplied with the Thin Cylinder apparatus. The guide includes full details of the equipment, detailed experiment procedures, theory and results.

For quick and reliable tests, TecEquipment can supply VDAS® which gives accurate real-time data capture, monitoring and display, calculation and charting of all important readings on a computer (computer not included).

Standard Features

- Supplied with comprehensive user guide
- Five-year warranty
- Manufactured in accordance with the latest European Union directives

Recommended Ancillaries

- Bench-mounted version of the Versatile Data Acquisition System (VDAS-B)

Experiments

Investigations into stresses and strains in a thin cylinder, to give students an understanding of:

- Longitudinal stress, hoop (or circumferential) stress, radial stress and biaxial stress
- The behaviour of the cylinder under both open and closed-end conditions
- The use of strain gauges
- The stress strain relationship and value of Young's modulus for the cylinder material
- Indirect strain and stress
- The value of Poisson's ratio for the cylinder material
- The use of Mohr's circle to calculate the shear strain at any position in the cylinder
- The use of the 'superposition method' to find the principal strains
- The effect of the biaxial stress system
- Sources of error in their experiments

Essential Services

Electrical supply:

Single-phase 100 VAC to 230 VAC, 50 Hz to 60 Hz

Operating Conditions

Operating environment: Laboratory

Storage temperature range: -25°C to +55°C (when packed for transport)

Operating temperature range: +5°C to +40°C

Operating relative humidity range: 80% at temperatures < 31°C decreasing linearly to 50% at 40°C

Sound Levels

The measured sound pressure level of this apparatus is less than 70 dB(A).

Specifications

Dimensions: 715 mm x 370 mm x 380 mm

Nett weight: 30 kg

Bench space needed: 1 m x 0.5 m

Oil: Shell Tellus S2 M32 or equivalent (supplied with 500 mL spare)

Cylinder and reservoir capacity: Approximately 2 litres

Cylinder material: Aluminium alloy

Pressure: Displayed by Bourdon gauge, with output for VDAS®

Strain: Measured by electrical resistance gauges and displayed digitally, with output for VDAS®

- TecEquipment Ltd, Bonsall Street, Long Eaton, Nottingham NG10 2AN, UK
- **T** +44 115 972 2611 • **F** +44 115 973 1520 • **E** info@tecequipment.com • **W** www.tecequipment.com
- An ISO 9001 certified company
- VDAS is a registered trademark of TecEquipment Ltd

tradition.

innovation.

integration.

infoWERK is a leading expert in the development of eLearning courseware, learning system solutions, teaching and AV equipment.

Furthermore infoWERK is the representative and system integrator of "TecQuipment".

TecQuipment is one of the global leaders in technical teaching equipment for engineering. If you are interested in one of TecQuipment's products feel free to contact us at:



infoWERK Medien & Technik GmbH

Martinsbühel 6 / A-6170 Zirl / Austria

Phone: +43 (0) 5238 52099-0 / Fax: +43 (0) 5238 52099-40

E-Mail: info@infowerk.at / Website: infowerk.at

Otto-Dürr-Straße 25

D-70435 Stuttgart, Zuffenhausen/ Germany

Phone: +49 (0) 711 342471-0 / Fax: +49 (0) 711 342471-11

E-Mail: info@de.infowerk.at / Website: infowerk.at