

TD1001

Gay-Lussac's Law

Shows the relationship between pressure and temperature of a fixed volume of ideal gas

Works with
VDAS[®]



Typical screenshot of the VDAS[®] software

- Self-contained bench-top module
- Demonstrates Gay-Lussac's Law relating pressure and temperature of an ideal gas (air)
- Simple and safe - needs no tools, uses low pressures and a thermally-insulated heater
- Includes thermocouples and a pressure sensor connected to a digital display
- Electronic controller to accurately regulate temperature
- Can connect to TecEquipment's Versatile Data Acquisition System (VDAS[®])

- TecEquipment Ltd, Bonsall Street, Long Eaton, Nottingham NG10 2AN, UK
- **T** +44 115 972 2611 • **F** +44 115 973 1520 • **E** info@tecquipment.com • **W** www.tecquipment.com
- An ISO 9001 certified company
- VDAS is a registered trademark of TecEquipment Ltd

TD1001

Gay-Lussac's Law

Description

The bench-mounting equipment includes a back plate that holds a low-pressure vessel. The vessel holds a fixed volume of air surrounded by an insulated heater, controlled by an electronic temperature controller.

A hand-operated valve at the bottom of the vessel allows students to normalize the air in the vessel to ambient conditions.

The equipment uses normal clean dry air, as it behaves as an ideal gas over the range of pressures used in this equipment.

A thermocouple measures the temperature of the heater surface for the controller. Two thermocouples measure the temperature of the air in the vessel. A pressure transducer measures the pressure of the heated air in the vessel. A digital display shows the absolute pressure, both temperatures and their average value.

Students set the controller for the range of temperatures needed during the experiment. They then record the changes in pressure as the temperature increases and plot the results to prove Gay-Lussac's Law.

The experiment can also work in reverse; students heat the vessel, open the valve to normalize the air in the vessel, then shut the valve. They then record the pressure and temperature drop as the vessel cools naturally. This gives a different starting point and results which will fall below local ambient. Due to the slow nature of natural cooling, the optional VDAS is helpful in this test to log results automatically.

You can do tests with or without a computer connected. However, for quicker tests with easier recording of results, TecEquipment can supply the optional Versatile Data Acquisition System (VDAS). This gives accurate real-time data capture, monitoring and display, calculation and charting of all the important readings on a computer (computer not included).

Standard Features

- Supplied with comprehensive user guide
- Five-year warranty
- Made in accordance with the latest European Union directives

Experiments

- Demonstrates change of pressure of a fixed volume of gas during heating.
- Proving Gay-Lussac's Law by experiment.
- The principle of a vapour pressure thermometer.

Recommended Ancillaries

- VDAS-B (bench-top version of the Versatile Data Acquisition System)

Operating Conditions

Operating environment:
Laboratory

Note: This equipment accurately measures temperatures and pressures with respect to normal atmospheric conditions. For best results, you must use it in a laboratory with a stable temperature around 20°C and away from direct heat sources.

Storage temperature range:
-25°C to +55°C (when packed for transport)

Operating temperature range:
+15°C to +30°C

Operating relative humidity range:
80% at temperatures < 31°C decreasing linearly to 50% at 40°C

Sound Levels

Less than 70 dB(A)

Essential Services

Bench Space Needed - 630 mm x 520 mm plus space for a suitable computer if you need to use the optional VDAS

Electrical Supply (determined by order):
110 VAC to 120 VAC or 220 VAC to 240 VAC
50 Hz to 60 Hz at 0.5 A

Technical Details

Nett Dimensions and Weight
630 mm x 520 mm x 600 mm high and 18 kg

Approximate Packed Volume and Weight
0.3 m³ and 20 kg

- TecEquipment Ltd, Bonsall Street, Long Eaton, Nottingham NG10 2AN, UK
- **T** +44 115 972 2611 • **F** +44 115 973 1520 • **E** info@tecequipment.com • **W** www.tecequipment.com
- An ISO 9001 certified company
- VDAS is a registered trademark of TecEquipment Ltd



tradition.

innovation.

integration.

infoWERK is a leading expert in the development of eLearning courseware, learning system solutions, teaching and AV equipment.

Furthermore infoWERK is the representative and system integrator of "TecQuipment".

TecQuipment is one of the global leaders in technical teaching equipment for engineering. If you are interested in one of TecQuipment's products feel free to contact us at:



infoWERK Medien & Technik GmbH

Martinsbühel 6 / A-6170 Zirl / Austria

Phone: +43 (0) 5238 52099-0 / Fax: +43 (0) 5238 52099-40

E-Mail: info@infowerk.at / Website: infowerk.at

Otto-Dürr-Straße 25

D-70435 Stuttgart, Zuffenhausen/ Germany

Phone: +49 (0) 711 342471-0 / Fax: +49 (0) 711 342471-11

E-Mail: info@de.infowerk.at / Website: infowerk.at