




Thermodynamics

TD360

Bench-top Heat Exchangers Service Module

Service module for TecEquipment's bench-top heat exchangers

Works with
VDAS[®]



Screenshot of the optional VDAS[®] software



- A bench-top service module with optional small-scale demonstration heat exchangers – designed for teaching
- Optional heat exchangers include the most common types used in industry (tubular, plate, shell and tube, and a jacketed vessel with coil and stirrer)
- Simple and safe to use – foolproof fittings allow students to change and connect the optional heat exchangers quickly and easily – needs no tools
- Clear digital displays of all readings – you do not need a computer to work it or take readings
- Can connect to TecEquipment's Versatile Data Acquisition System (VDAS[®])


- TecEquipment Ltd, Bonsall Street, Long Eaton, Nottingham NG10 2AN, UK
- **T** +44 115 972 2611 • **F** +44 115 973 1520 • **E** info@tecquipment.com • **W** www.tecquipment.com
- An ISO 9001 certified company
- VDAS is a registered trademark of TecEquipment Ltd



TD360

Bench-top Heat Exchangers Service Module

Description



Shown fitted with one of the optional heat exchanger modules and optional VDAS-F Interface

The Service Module (TD360) is the core of the TD360 range. It provides hot and cold water to the heat exchangers and all the instruments needed to measure their performance. All fluid connections to the optional heat exchangers are self-sealing quick connectors – for safety and simplicity. The hot and cold fluid streams have different connectors to reduce errors.

The services module's hot water system includes a tank with a PID controlled electric heater, a pump and tank level indicators. An electrically operated valve opens to let water in to fill the tank. The tank has protection in case of over temperature, low water level and overfilling. The hot water system gives stable flow rates and temperatures.

The services module's cold water circuit has a flow regulator and connection for an external mains water supply.

Both the cold and hot water system have precision needle valves and turbine flow meters to control and measure the flow rates.

Thermocouples at the connectors measure hot and cold inlet and outlet fluid stream temperatures. Some heat-exchangers also have built-in thermocouples for extra temperature measurements. Clear, multiline digital displays show the temperatures and flow rates of the fluid streams.

All optional heat exchangers have the same nominal heat transfer area and wall thickness, so students can compare them directly.

You can do tests with or without a computer connected. However, for quicker tests with easier recording of results, TecQuipment can supply the optional Versatile Data Acquisition System (VDAS®). This gives accurate real-time data capture, monitoring and display, calculation and charting of all the important readings on a computer (computer not included).

Available Experiment Modules

- Concentric Tube Heat Exchanger (TD360a)
- Plate Heat Exchanger (TD360b)
- Shell and Tube Heat Exchanger (TD360c)
- Jacketed Vessel with Coil and Stirrer (TD360d)

Note: You need at least one of the optional heat exchangers to do experiments. TecQuipment recommends that you buy the Concentric Tube Heat Exchanger (TD360a) first, because it has extra temperature measuring points.

Recommended Ancillaries

- VDAS-F (frame-mounted version of the Versatile Data Acquisition System)

Standard Features

- Supplied with comprehensive user guide
- Five-year warranty
- Made in accordance with the latest European Union directives

Operating Conditions

Operating environment:
Laboratory

Storage temperature range:
–25°C to +55°C (when packed for transport)

Operating temperature range:
+5°C to +40°C

Operating relative humidity range:
80% at temperatures < 31°C decreasing linearly to 50% at 40°C

Sound Levels

Less than 70 dB(A)

Essential Services

Bench space needed: 800 mm x 700 mm

Electrical supply: Single-phase 220 VAC 50 Hz or two-phase 220 VAC 60 Hz at 13 A

Clean water supply and waste: 5 L.min⁻¹ at a minimum 1 bar and maximum 3 bar.

Specifications

Dimensions and weight:
760 mm wide x 610 mm front to back x 600 mm high and 40 kg

- TecQuipment Ltd, Bonsall Street, Long Eaton, Nottingham NG10 2AN, UK
- **T** +44 115 972 2611 • **F** +44 115 973 1520 • **E** info@tecquipment.com • **W** www.tecquipment.com
- An ISO 9001 certified company
- VDAS is a registered trademark of TecQuipment Ltd

tradition.

innovation.

integration.

infoWERK is a leading expert in the development of eLearning courseware, learning system solutions, teaching and AV equipment.

Furthermore infoWERK is the representative and system integrator of "TecQuipment".

TecQuipment is one of the global leaders in technical teaching equipment for engineering. If you are interested in one of TecQuipment's products feel free to contact us at:



infoWERK Medien & Technik GmbH

Martinsbühel 6 / A-6170 Zirl / Austria

Phone: +43 (0) 5238 52099-0 / Fax: +43 (0) 5238 52099-40

E-Mail: info@infowerk.at / Website: infowerk.at

Otto-Dürr-Straße 25

D-70435 Stuttgart, Zuffenhausen/ Germany

Phone: +49 (0) 711 342471-0 / Fax: +49 (0) 711 342471-11

E-Mail: info@de.infowerk.at / Website: infowerk.at