


TE19

Thermal Conductivity Experiment

Shows heat conduction in different metals at ambient temperature

Works with
VDAS®



Screenshot of the optional VDAS® software

- Bench-mounting apparatus that demonstrates heat conduction along metal test specimens
- Uses low temperatures to reduce heat losses, for accurate results
- Test area pressure can be reduced to a vacuum for accurate results (needs optional vacuum pump)
- Supplied with a selection of different metal specimens
- Specimens are reusable and easily changed
- Connects to TecQuipment's optional Versatile Data Acquisition System (VDAS®) for automatic data acquisition

- TecQuipment Ltd, Bonsall Street, Long Eaton, Nottingham NG10 2AN, UK
- **T** +44 115 972 2611 • **F** +44 115 973 1520 • **E** info@tecquipment.com • **W** www.tecquipment.com
- An ISO 9001 certified company
- VDAS is a registered trademark of TecQuipment Ltd

TE19

Thermal Conductivity Experiment

Description

The Thermal Conductivity Experiment (TE19) is in two parts:

1) A base unit with a Peltier cooler, a vacuum gauge and a moulded glass dome (vacuum vessel).

2) A separate control module, with heater and cooler controls and a multiline display. It also includes a connection for TecEquipment's optional VDAS® and a suitable computer (not supplied) for automatic data acquisition.

Included with the equipment are specimens of different metals. Each specimen has a small resistive heater at one end, the other end clamps to the cooler of the base unit. The heater at one end and the cooler at the other give a controlled heat flow along the specimen. Attached to the specimen are two thermocouples at a precise distance apart, near to each end. Students use the temperature difference between the thermocouples to find the thermal conductivity of the material. They then compare it with other materials and with given values.

The glass dome covers the area around the specimen. The students may use a suitable vacuum pump (not supplied) to remove the air inside the dome. This reduces heat loss due to convection, to give more accurate results.

For quick and reliable tests, TecEquipment can supply the optional VDAS® (Versatile Data Acquisition System). VDAS® gives accurate real-time data capture, monitoring and display, calculation and charting of all important readings on a computer (computer not included).

Standard Features

- Supplied with comprehensive user guide
- Five-year warranty
- Manufactured in accordance with the latest European Union directives

Essential Ancillaries

- Laboratory Vacuum Pump (RE19)

Recommended Ancillaries

- VDAS-B

A bench-mounted version of the Versatile Data Acquisition System

Experiments

- Thermal conductivity of different metals

Essential Services

Electrical supply:

Single-phase, 230 V/110 V, 50/60 Hz, 0.5 A

Bench space needed:

Solid, level bench top 1300 mm x 400 mm, plus space for the optional vacuum pump.

Operating Conditions

Operating environment:

Laboratory

Storage temperature range:

-25°C to +55°C (when packed for transport)

Operating temperature range:

+5°C to +40°C

Operating relative humidity range:

80% at temperatures < 31°C decreasing linearly to 50% at 40°C

Sound Levels

TE19 – Less than 70 dB(A)

Refer to the manufacturer's details for sound levels of any vacuum pump that you use with this equipment.

Specifications

Nett dimensions and weight:

Base Unit: 650 mm x 350 mm x 600 mm and 20 kg

Control Module: 600 mm x 350 mm x 180 mm and 11 kg

Heater:

Electrical resistance, fitted at the end of each specimen

Cooler:

Peltier device with heat sink

Thermocouples:

K-type

Vacuum gauge:

0 to -1 bar

Specimen metals:

- Brass
- Aluminium
- Copper

- TecEquipment Ltd, Bonsall Street, Long Eaton, Nottingham NG10 2AN, UK
- **T** +44 115 972 2611 • **F** +44 115 973 1520 • **E** info@tecequipment.com • **W** www.tecequipment.com
- An ISO 9001 certified company
- VDAS is a registered trademark of TecEquipment Ltd



TecEquipment

tradition.

innovation.

integration.

infoWERK is a leading expert in the development of eLearning courseware, learning system solutions, teaching and AV equipment.

Furthermore infoWERK is the representative and system integrator of "TecQuipment".

TecQuipment is one of the global leaders in technical teaching equipment for engineering. If you are interested in one of TecQuipment's products feel free to contact us at:



infoWERK Medien & Technik GmbH

Martinsbühel 6 / A-6170 Zirl / Austria

Phone: +43 (0) 5238 52099-0 / Fax: +43 (0) 5238 52099-40

E-Mail: info@infowerk.at / Website: infowerk.at

Otto-Dürr-Straße 25

D-70435 Stuttgart, Zuffenhausen/ Germany

Phone: +49 (0) 711 342471-0 / Fax: +49 (0) 711 342471-11

E-Mail: info@de.infowerk.at / Website: infowerk.at