



Control Engineering

TE3300/05

Temperature Process Training System

For a wide range of practical experiments in temperature control



- Shows automatic control of temperature using proportional, proportional plus integral, and proportional, integral plus derivative (PID) control
- Uses industry-standard parts to make it ideal for industrial, vocational and academic training
- Shows operation, calibration and tuning of temperature transmitters and thermocouples
- Compact, mobile and fully self-contained
- Includes delay coil to mimic realistic time lag due to a process
- Connects to the TE3300/06 Computer Control System for distributed control
- Safe, practical and realistic

- TecEquipment Ltd, Bonsall Street, Long Eaton, Nottingham NG10 2AN, UK
- **T** +44 115 972 2611 • **F** +44 115 973 1520 • **E** info@tecequipment.com • **W** www.tecequipment.com
- An ISO 9001 certified company

TE3300/05

Temperature Process Training System

Description

The Temperature Process Training System is a compact and mobile unit for a wide range of experiments in temperature control. It gives students a greater understanding of the stability of simple control systems.

The self-contained unit can do many experiments, but it can also connect to the optional Computer Control System (TE3300/06) for distributed control.

The main parts of the Temperature Process Training system are:

- Industrial controller with autotune feature
- Two-channel chart recorder
- Heat-exchanger and fan
- Temperature transmitter
- Thermocouples
- Delay coil
- In-line heater
- Three-speed pump
- Reservoir

To perform experiments, students fill the reservoir with clean water and prime the system. They then set the controller to regulate the power to the in-line heater and control the temperature of the water at any of three places. The heat-exchanger removes the heat from the water, to give quicker experiments. The thermocouples (selected by a three-way switch) give feedback to the controller.

For a realistic experience, the equipment has industrial-standard instrumentation and parts.

The apparatus includes one gate valve that works as a flow bypass. A chart recorder shows and logs the changes of the process variable (temperature) and the controller output.

Note: The chart recorder is paperless, so you need a suitable computer and colour printer if you need to print out hard copies of the chart recorder traces.

A socket on the side of the apparatus links to the Computer Control System (TE3300/06, available separately).

Standard Features

- Supplied with comprehensive user guides
- Five-year warranty
- Made in accordance with the latest European Union directives

Recommended Ancillaries

- Stopwatch (SW1)
- Thermocouple calibrator*
- One-litre container with accurate scale and accurate weighing machine (to measure water flow rate)*
- Computer Control System (TE3300/06)

* not supplied by TecQuipment

Experiments

- Proportional, integral and derivative control
- Setting up and demonstrating automatic control
- The principles of loop control and the calibration and tuning of temperature transmitters and thermocouples
- Operation of a temperature control system
- Distributed control (when used with the TE3300/06 Computer Control System)

Operating Conditions

Operating environment:
Laboratory environment

Storage temperature range:
-25°C to +55°C (when packed for transport)

Operating temperature range:
+5°C to +40°C

Operating relative humidity range:
80% at temperatures < 31°C decreasing linearly to 50% at 40°C

Sound Levels

Less than 70 dB (A)

Essential Services

Electrical supply:
Single-phase 230 V at 3 A or 110 V at 7 A, 50/60 Hz (determined at order)

Dimensions and Weights

Nett: 700 mm x 800 mm x 1750 and 120 kg

Packed: Approximately 1 m³ and 150 kg

- TecQuipment Ltd, Bonsall Street, Long Eaton, Nottingham NG10 2AN, UK
- **T** +44 115 972 2611 • **F** +44 115 973 1520 • **E** info@tecquipment.com • **W** www.tecquipment.com
- An ISO 9001 certified company

tradition.

innovation.

integration.

infoWERK is a leading expert in the development of eLearning courseware, learning system solutions, teaching and AV equipment.

Furthermore infoWERK is the representative and system integrator of "TecQuipment".

TecQuipment is one of the global leaders in technical teaching equipment for engineering. If you are interested in one of TecQuipment's products feel free to contact us at:



infoWERK Medien & Technik GmbH

Martinsbühel 6 / A-6170 Zirl / Austria

Phone: +43 (0) 5238 52099-0 / Fax: +43 (0) 5238 52099-40

E-Mail: info@infowerk.at / Website: infowerk.at

Otto-Dürr-Straße 25

D-70435 Stuttgart, Zuffenhausen/ Germany

Phone: +49 (0) 711 342471-0 / Fax: +49 (0) 711 342471-11

E-Mail: info@de.infowerk.at / Website: infowerk.at